

# **Introduction to web-based bioinformatics tools - Resfinder, Virulence Finder using CGE**

**Dr Mabel Aworh-Ajumobi**  
**North Carolina State University, Raleigh**


# Learning objectives:

- To learn how to access CGE Tools
- To use ResFinder in determination of AMR genes
- To use Virulence Finder in identification of virulence genes in bacteria species
- To retrieve and interpret outputs from ResFinder and Virulence Finder

# Bioinformatics web-based tools

- Variety of methods and tools are available to analyze bacterial pathogens
- Most bioinformatics tools are implemented in Unix environments
- Require at least some bioinformatics expertise for usage
- Web-based bioinformatics tools
  - Often free for use
  - Do not require computational power from the user
  - Limited bioinformatics knowledge
  - In some cases demand that that users deposit the analyzed data in public repositories
  - Always make an effort to browse through the documentation of web based platforms
    - Helps with choice of parameters and interpretation of results

# Examples of free online web- based tools



# Accessing the CGE homepage


## Center for Genomic Epidemiology

Home

Services

Publications

Contact



### News

ResFinder – an open online resource for identification of antimicrobial resistance genes in next-generation sequencing data and prediction of phenotypes from genotypes

January 2022

[Link to article....](#)

MINTyper: an outbreak-detection method for accurate and rapid SNP typing of clonal clusters with noisy long reads

April 2021

[Link to article....](#)

Automated download and clean-up of family specific databases for kmer-based virus identification

October 2020

[Link to article....](#)

CRHP Finder, a webtool for the detection of

# Overview of CGE Tools

## Phenotypic

### [ResFinder](#)

Identification of acquired antibiotic resistance genes.

### [ResFinderFG](#)

Identification of functional metagenomic antibiotic resistance determinants.

### [LRE-finder](#)

Identification of genes and mutations leading to linezolid resistance.

### [KmerResistance](#)

Identification of acquired antibiotic resistance genes using Kmers.

### [PathogenFinder](#)

Prediction of a bacteria's pathogenicity towards human hosts.

### [VirulenceFinder](#)

Identification of acquired virulence genes.

### [Restriction-ModificationFinder](#)

Determination of Restriction-Modification sites (based on [REBASE](#).)

### [SPIFinder](#)

SPIFinder identifies Salmonella Pathogenicity Islands.

### [ToxFinder](#)

ToxFinder identifies genes involved in mycotoxin synthesis.

## Phylogeny

### [MINTyper](#)

Identification of SNPs with automatic filtering, masking and site validation together with inferred phylogeny based on both long and short sequencing data.

### [CSIPhylogeny](#)

CSI Phylogeny calls SNPs, filters the SNPs, does site validation and infers a phylogeny based on the concatenated alignment of the high quality\* SNPs.

### [NDtree](#)

NDtree constructs phylogenetic trees from Single-End or Pair-End FASTQ files.

### [Evergreen](#)

Evergreen generates a forest of constantly updated phylogenetic trees with publicly available whole-genome sequencing data from foodborne, bacterial isolates that were deposited in the short sequencing read archives (NCBI SRA/ENA).

### [TreeViewer](#)

Phylogeny Tree Viewer.

## Metagenomics

### [CCMetagen](#)

CCMetagen: Comprehensive and accurate identification of eukaryotes and prokaryotes in metagenomic data.

# ResFinder – for detecting resistance genes

- Is based on curated database, public databases as well as scientific papers
- Is web-based friendly interface and freely accessible tool
- Will detect the presence of whole genome resistance genes
- The database contains genes for 14 major antimicrobial drug classes

We want to learn  
how to detect AMR  
genes in bacteria  
genome

## ResFinder 4.1

ResFinder identifies acquired genes and/or finds chromosomal mutations mediating antimicrobial resistance in total or partial DNA sequence of bacteria.

The database is curated by:  
**Frank Møller Aarestrup**  
([click to contact](#))

### Updates

ResFinder and PointFinder software: ([2020-10-21](#))

ResFinder database: ([2021-01-20](#))

PointFinder database: ([2019-07-02](#))

Chromosomal point mutations ☐

Acquired antimicrobial resistance genes ☐

### Select species


Escherichia coli\*

\*Chromosomal point mutation database exists

### Select type of your reads

Assembled Genome/Contigs

If you get an "Access forbidden. Error 403": Make sure the start of the web address is https and not just http. Fix it by clicking [here](#).

 Isolate File

Name	Size	Progress	Status

 Upload

 Remove



# To upload your sequence file

## ResFinder 4.1


ResFinder identifies acquired genes and/or finds chromosomal mutations mediating antimicrobial resistance in total or partial DNA sequence of bacteria.

The database is curated by:  
**Frank Møller Aarestrup**  
(click to contact)


### Updates

ResFinder and PointFinder software: (2020-10-21)  
ResFinder database: (2021-01-20)  
PointFinder database: (2019-07-02)

---

Chromosomal point mutations ☐ 


---

Acquired antimicrobial resistance genes ☐ 

---

Select species  


Escherichia coli\*



\*Chromosomal point mutation database exists

Select type of your reads  

Assembled Genome/Contigs



If you get an "Access forbidden. Error 403": Make sure the start of the web address is https and not just http. Fix it by clicking [here](#).

Isolate File

Name	Size	Progress	Status

Upload

Remove

Ensure to check all the important boxes has shown by the green tick marks

# Assessing your results

## Center for Genomic Epidemiology

### Your job is being processed

Wait here to watch the progress of your job, or fill in the form below to get an email message upon completion.

To get notified by email:

This page will update itself automatically.



● cgewebface@cbs.dtu.dk

To: maworh@gmail.com



Fri, Feb 3 at 5:55 PM

Dear sir/madam,

Your job 63DDB57D0000051B70E4CBE5 has finished on our ResFinder-4.1 server.

You may visit the following link to study the results:

<https://cge.food.dtu.dk/cgi-bin/webface.fcgi?jobid=63DDB57D0000051B70E4CBE5>

The data will expire at 2023-02-11 02:31:41.

Yours sincerely,  
CGE webserver

You can wait for the page to update after analysis is completed or choose to get notified by email

# ResFinder-4.1 Server - Results

Input Files: *isolate2.fasta*

**Warning:**

One or more resistance genes does not exist in the phenotype database. The Summary table does not take this into account.

escherichia coli complete			
Antimicrobial	Class	WGS-predicted phenotype	Genetic background
amikacin	aminoglycoside	No resistance	blaTEM-1B (blaTEM-1B_AY458016)
tigecycline	tetracycline	No resistance	
tobramycin	aminoglycoside	No resistance	
cefepime	beta-lactam	No resistance	
chloramphenicol	amphenicol	No resistance	
piperacillin+tazobactam	beta-lactam	No resistance	
cefoxitin	beta-lactam	No resistance	
ampicillin	beta-lactam	Resistant	
ampicillin+clavulanic acid	beta-lactam	No resistance	
cefotaxime	beta-lactam	No resistance	
ciprofloxacin	quinolone	Resistant	gyrA (p.S83L)
colistin	polymyxin	No resistance	
sulfamethoxazole	folate pathway antagonist	Resistant	sul2 (sul2_AY034138)
imipenem	beta-lactam	No resistance	dfrA14 (dfrA14_KF921535)
trimethoprim	folate pathway antagonist	Resistant	
nalidixic acid	quinolone	Resistant	
ertapenem	beta-lactam	No resistance	gyrA (p.S83L)
tetracycline	tetracycline	Resistant	
fosfomycin	fosfomycin	No resistance	
tetracycline	tetracycline	No resistance	tet(A) (tet(A)_AJ517790)
fosfomycin	fosfomycin	No resistance	

Aminoglycoside									
Resistance gene	Identity	Alignment Length/Gene Length	Position in reference	Contig or Depth	Position in contig	Phenotype	PMID	Accession no.	Notes
aph(6)-Id	100.0	837/837	1..837	NODE_79_length_7882_cov_4.782979	3194..4030	streptomycin	2653965	<a href="#">M28829</a>	Alternative name strB
aph(3'')-Ib	100.0	804/804						<a href="#">AF321551</a>	Alternative name strA

**Alignment length** is the length between the best matching resistance gene and corresponding sequence in the genome while **Gene length** is the length of the best matching resistance gene in the database

Resistance gene	Identity	Alignment Length/Gene Length	Position in reference	Contig or Depth	Position in contig	Phenotype	PMID	Accession no.	Notes
tet(A)	100.0	1200/1200	1..1200	NODE_79_length_7882_cov_4.782979	387..1586	doxycycline,tetracycline	12654659	<a href="#">AJ517790</a>	

The dark green colour indicates a perfect match for a given gene. The % ID is 100 and means that the sequence in the genome covers the entire length of the resistance gene in the database

Beta-lactamase									
Resistance gene	Identity	Alignment Length/Gene Length	Position in reference	Contig or Depth	Position in contig	Phenotype	PMID	Accession no.	Notes
blaTEM-1B	100.0	861/861	1..861	NODE_102_length_2700_cov_6.467548	1323..2183	amoxicillin,ampicillin,cephalothin,piperacillin,ticarcillin	15388431	<a href="#">AY458016</a>	Class A

# Other ways to access your results

## Download acquired AMR gene results:

[Results as text](#)[Hit in genome sequences](#)[Resistance gene sequences](#)[Results as tabseperated file](#)

## Download Chromosomal point mutation results:

[Results as tabseperated file](#)[Results as a text file](#)

**Selected %ID threshold for ResFinder: 90 %**

**Selected minimum length for ResFinder: 60 %**

**Selected %ID threshold for PointFinder: 90 %**

**Selected minimum length for PointFinder: 60 %**

[extended output](#)

# Click on extended output for more information

extended output

dfrA14, ID: 100.00 %, Alignment Length/Gene Length: 474/474, Coverage: 100.0, Positions in reference: 1..474, Contig name: NODE\_104\_length\_2160\_cov\_

Resistance gene seq: TTGAAAGTATCATTGATGGCTGCGAAAAGCGAAAAACGGCGTGATTGGTTGCGGTCCAGAC  
Hit in genome: TTGAAAGTATCATTGATGGCTGCGAAAAGCGAAAAACGGCGTGATTGGTTGCGGTCCAGAC

Resistance gene seq: ATACCCTGGTCCGCGAAAGGGGAGCAGCTACTTTTTAAAGCATTGACCTACAATCAGTGG  
Hit in genome: ATACCCTGGTCCGCGAAAGGGGAGCAGCTACTTTTTAAAGCATTGACCTACAATCAGTGG

Resistance gene seq: CTTCTGGTGGGTCGCAAGACGTTTGAATCTATGGGCGCACTCCCAATAGGAAATACGCG  
Hit in genome: CTTCTGGTGGGTCGCAAGACGTTTGAATCTATGGGCGCACTCCCAATAGGAAATACGCG

Resistance gene seq: GTCGTTACCCGCTCAGGTTGGACATCAAATGATGACAATGTAGTTGTATTTTCAGTCAATC  
Hit in genome: GTCGTTACCCGCTCAGGTTGGACATCAAATGATGACAATGTAGTTGTATTTTCAGTCAATC

Resistance gene seq: GAAGAGGCCATGGACAGGCTAGCTGAATTCACCGGTCACGTTATAGTGTCTGGTGGCGGA  
Hit in genome: GAAGAGGCCATGGACAGGCTAGCTGAATTCACCGGTCACGTTATAGTGTCTGGTGGCGGA

Resistance gene seq: GAAATTTACCGAGAAACATTACCCATGGCCTCTACGCTCCACTTATCGACGATCGACATC  
Hit in genome: GAAATTTACCGAGAAACATTACCCATGGCCTCTACGCTCCACTTATCGACGATCGACATC

Resistance gene seq: GAGCCAGAGGGGGATGTTTTCTTCCCGAGTATTCCAAATACCTTCAAGTTGTTTTTGAG

This shows you the alignment between your sequence in the genome and the resistance genes. Green colour shows a perfect match between the sequence in the genome and resistance genes in the database

# To detect virulence genes

## VirulenceFinder 2.0

[Service](#) [Instructions](#) [Output](#) [Article abstract](#) [Citations](#) [Version history](#)

Software version: 2.0.3 (2020-05-21)  
Database version: (2022-12-02)

The database is curated by:  
**Flemming Scheutz, SSI**  
(click to contact)

**Select species**

Listeria

S. aureus

Escherichia coli

Enterococcus

**Select threshold for %ID**

90 %

**Select type of your reads**

Only data from one single isolate should be uploaded. If raw sequencing reads are uploaded KMA will be used for mapping. KMA supports the following sequencing platforms: Illumina, Ion Torrent, Roche 454, SOLiD, Oxford Nanopore, and PacBio.

Assembled or Draft Genome/Contigs\* ( v )



Choose File(s)

Name	Size	Progress	Status

Upload Remove

This tool is used to analyse bacterial sequence data to detect the presence of virulence genes





Choose file(s)

Name	Size	Progress	Status

Upload Remove



# Assessing your results

## Center for Genomic Epidemiology

### Your job is being processed

Wait here to watch the progress of your job, or fill in the form below to get an email message upon completion.  
**maworh@gmail.com**

To get notified by email:

This page will update itself automatically



● **CGEwebface@cbs.dtu.dk** <cgewebface@cbs.dtu.dk>  
To: maworh@gmail.com



Fri, Feb 3 at 5:37 PM

Dear sir/madam,

Your job 63DDB14700007D13AF2F51DF has finished on our VirulenceFinder-2.0 server.

You may visit the following link to study the results:

<https://cge.food.dtu.dk/cgi-bin/webface.fcgi?jobid=63DDB14700007D13AF2F51DF>

The data will expire at 2023-02-11 02:13:43.

Yours sincerely,  
CGE webserver

# VirulenceFinder-2.0 Server - Results

Organism(s): *Escherichia coli*

The grey color indicates a warning due to a non-perfect match; Query length is shorter than template length

Virulence genes for Escherichia coli						
Virulence factor	Identity	Query / Template length	Contig	Position in contig	Protein function	Accession number
AsIA	93.9	1657 / 1656	NODE_13_length_88665_cov_5.581852	37650..39305		<a href="#">CP054236</a>
csgA	100	456 / 456	NODE_89_length_9640_cov_5.145275	294..749	curlin major subunit CsgA	<a href="#">CP069646</a>
fdeC	92.88	4258 / 4254	NODE_14_length_87312_cov_5.365935	1478		
fimH	99.59	489 / 489	NODE_68_length_19341_cov_4.739565	6213		
hlyE	100	918 / 918	NODE_5_length_144015_cov_5.105777	104334..105251	Avian E. coli haemolysin	<a href="#">ECU57430</a>
nlpl	99.77	885 / 885	NODE_12_length_94487_cov_5.441109	41575..42459	lipoprotein Nlpl precursor	<a href="#">CP000243</a>

The dark green color indicates a perfect match for a given gene. The % ID is 100 and means that the sequence in the genome covers the entire length of the virulence gene in the database

Shiga-toxin genes						
Virulence factor	Identity	Query / Template length	Contig	Position in contig	Protein function	Accession number
No hit found						

The light green color indicates a warning due a non-perfect match, ID < 100% although both query and template lengths are equal

The red color indicates that no virulence gene was detected for the shiga toxin

# Click on extended output for more information

extended output

Results as text

Results tsv

Hits in genome seqs

Virulence factor seqs

**Input Files:** *ma\_011.fasta*

## CITATIONS

For publication of results, please cite:

- Real-time whole-genome sequencing for routine typing, surveillance, and outbreak detection of verotoxigenic Escherichia coli. Joensen KG, Scheutz F, Lund O, Hasman H, Kaas RS, Nielsen EM, Aarestrup FM. J. Clin. Microbiol. 2014. 52(5): 1501-1510.  
View the [abstract](#)

## *fimH\_NA*

template	TTCGCCTGTAAAACCGCCAATGGTACAGCTATCCCTATTGGCGGTGGCAGCGCTAATGTT
query	TTCGCCTGTAAAACCGCCAATGGTACAGCTATCCCTATTGGCGGTGGCAGCGCTAATGTT
template	TATGTAAACCTTGCGCCCGCCGTGAATGTGGGGCAAACCTGGTCGTAGATCTTTCGACG
query	TATGTAAACCTTTTCGCCCGCCGTGAATGTGGGGCAAACCTGGTCGTAGATCTTTCGACG
template	CAAATCTTTTGCCATAACGATTATCCGGAAACCATTACAGACTATGTCACACTGCAACGA
query	CAAATCTTTTGCCATAACGATTATCCGGAAACCATTGACAGACTATGTCACACTGCAACGA
template	GGCTCGGCTTATGGCGGCGTGTTATCTAATTTTTCCGGGACCGTAAAATATAGTGGCAGT
query	GGCTCGGCTTATGGCGGCGTGTTATCTAATTTTTCCGGGACCGTAAAATATAGTGGCAGT
template	AGCTATCCATTTCTACCACCAGCGAAACGCCGCGCGTTGTTTATAATTCGAGAACGGAT
query	AGCTATCCATTTCTACCACCAGCGAAACGCCGCGCGTTGTTTATAATTCGAGAACGGAT
template	AAGCCGTGGCCGGTGGCGCTTTATTTGACGCCTGTGAGCAGTGCGGGCGGGGTGGCGATT
query	AAGCCGTGGCCGGTGGCGCTTTATTTGACGCCTGTGAGCAGTGCGGGCGGGGTGGCGATT
template	AAAGCTGGCTCATTAAATTGCCGTGCTTATTTTGCAGACAGACCAACAACCTATAACAGCGAT
query	AAAGCTGGCTCATTAAATTGCCGTGCTTATTTTGCAGACAGACCAACAACCTATAACAGCGAT
template	GATTTCCAGTTTGTGTGGAATATTTACGCCAATAATGATGTGGTGGTGCCCACTGGCGGC
query	GATTTCCAGTTTGTGTGGAATATTTACGCCAATAATGATGTGGTGGTGCCCACTGGCGGC
template	TGTGATGTT
query	TGTGATGTT

This shows you the alignment between your sequence in the genome & the virulence gene. The red colour shows a variation has occurred e.g instead of a base “G” there is a “T” and instead of base “A” we have a “G”

## *hlyE\_ECU57430*

# Now you should be able:

- To access ResFinder & VirulenceFinder and upload your sequences
- To access your results or outputs from ResFinder & VirulenceFinder
- To interpret the outputs from ResFinder & VirulenceFinder



*Thank you!*

Happy to take any  
questions.

You can contact me:  
**[maworh@gmail.com](mailto:maworh@gmail.com)**